

Sepamatic SL

Precise separation of blood components



Innovation Sepamatic SL

- Adaptive sensors
- Three-sided access to bag hanging
- Removal of air bubbles from plasma bag



Display with user friendly menu navigation



Innovation:
Pressing of air out
of the blood bag



Well designed slots for
hanging of the tubes



Sepamatic SL



Easy three-sided access to bag hanging



Sepamatic SL

- Triple press system
- Clean cellfree separation
- LAN Data transfer
- Open software with multiple functions
- Flow control system
- Softfilter press system

Features

- Easy handling
- Easy installation of new blood bag systems
- Adaptive sensors in the sealing clamps
- LAN - Ethernet data transfer
- WLAN data transfer
- Balances for whole blood, plasma, red cells and buffy coat
- Donation ID-System
- Individual application of different blood bags



Lmb

Innovation: Adaptive sensors for each tube and bag system



Address

Möslstr. 17
D-85445 Schwaig

Phone

+49 - 8122 - 880 96 - 0

Fax

+49 - 8122 - 880 96 - 60

Internet

www.lmb.de

E-mail

info@lmb.de

Sepamatic SL

Technical Specifications

■ Power Supply	220 VAC 10 %
■ Frequency	50 Hz
■ Power Consumption	300 VA
■ Fuse	2 x 4A / 250V
■ Dimensions	480 x 480 x 520 mm
■ Weight	35 kg
■ Balances	4
■ Clearance of balances	10 kg
■ Sealing clamps	5
■ Sealing power	90W / 50 Ohm
■ Sealing time	0,5 - 5 Sec.
■ Flow regulation	first device, from 0 to max
■ Press A - Movement clearance	100 mm
■ Press B - Two positions	2 x 50 mm
■ Data transfer	Wireless LAN
■ Sensor system	
■ Reading performance red cells in mother bag	6 Sensors in a row
■ Red cell sensor (Exit mother bagl)	4 pieces

Features

- Use of different blood bag systems
- Up to 5 bags are treated in one process
- 4 Balances
- Triple press system
- Top & Top and Top & Bottom Process
- Plasma filtering by special soffilter

Features

- Easy handling
- All applications possible
- Barcode reader
- Smart label technology
- Wireless LAN data transfer
- Easy connection to local PC network
- Data collection and data management software



Lmb

Address

Möslstr. 17
D-85445 Schwaig

Phone

+49 - 8122 - 880 96 - 0

Fax

+49 - 8122 - 880 96 - 60

Internet

www.lmb.de

E-mail

info@lmb.de